## **DATA SHEET**

PRODUCT CODE **DS41** 

## PREMIUM DOUBLE-SIDED TISSUE TAPE

Hand tearable double-sided tissue tape with good initial grab to a variety of surfaces. Highly conformable, with excellent temperature and solvent resistance.

## **FEATURES**

The product is imported on a printed liner: 'Remove to Expose Adhesive' – making it suitable for print finishing applications.

The product is available in a variety of widths, in both 50 meter hand rolls and 250 meter jumbo rolls.

PROPERTIES	UNIT	VALUE	TEST METHOD
Thickness	mm	0.11 (+/- 0.015)	ASTM D-1000
Adhesive		Acrylic	ASTM D-1000
Peel Adhesion to Steel	Kg/25mm	>0.60	ASTM D-1000
Holding Power	min	>700	PSTC-7
Tensile Strength	Kg/25mm	>1.5	ASTM D-1000
Storage Life	Months (temp)	12 at 23°C	ASTM D-1000

## **Terms & Conditions**

This product should be applied to a clean and dry surface. The adhesive of the product must not be exposed direct to the sunlight when applied. The above figures represent average values. The selection of a tape is the responsibility of the buyer. Tape should be stored at room temperature and should not be exposed to sunlight. We reserve the right to make changes without notice to our product and specifications

PH: 03 9544 1977 FAX: 03 9544 7113 sales@premiertapes.com.au **premiertapes.com.au** 



